

## BodyTone Lab Horizontal

The latest generation of body shaping device



LASERVILLE Medical Equipment Trading LLC

## 




#### BodyTone Lab Horizontal

#### Combining cardio workout with skin care treatment

The BodyTone Lab Horizontal is a unique device for weight loss and cellulite reduction. It combines physical exercise on a recumbent bike with a vacuum massage technology, infrared heating and collagen light applications, which in turn help in weight loss, burning fat, improving blood flow, maintaining a good health and rejuvenating the skin.

BODYT



#### **Treatments** Beyond fitness

Increasing metabolism

> Reducing cellulites

Recovering muscles

Avoiding spine & back pain

Burning fats

Reducing joint pain

#### There is much more to it!



#### **Technology** The evolution of cardio bike

BodyTone Lab Horizontal combines physical exercise on a recumbent bike with a vacuum massage technology, infrared heating and collagen.

## INFRARED

When the skin is exposed to the Infrared, it penetrates the fat cells and temporarily perforates their membranes, causing the lipids to leak out, in addition to the release of glycerol and water; i.e. breaking down fat more effectively and faster.

# 02

### VACUUM

Once the fat cells are broken down and the fatty acids, glycerol and water are released, the Vacuum will help the blood to flow into the parts of stubborn fat to use the fat for energy, and it will drain the water with the fatty acid and glycerol in the process.

# 03

## COLLAGEN

The Collagen simulates the body to reproduce collagen to help the skin to tighten, have more elasticity and smoothness and get rid of any sagginess.

#### **Benefits**

#### for your business



#### Functionality

BodyTone Lab Horizontal can be set up and prepared easily and quickly by any technician or personal trainer.



#### Design

It has a modern and luxurious design with a futuristic look and feel to make your studio or center look unique, stylish and appealing to your client.

#### Innovation

With BodyTone Lab Horizontal advanced technology and features, you can set up the parameters and customise the program based on your client's needs.



#### Training & Service

We offer assembly, installation and full service as well as, training.



#### ROI

Delivers high ROI and increases revenue with the increase of patient satisfaction

## for your client



#### Comfortability

The patient can enjoy a luxurious exercise experience on a relaxed bike chair in a comfortable lying down position



#### Less Time

In only 30 minutes session, your client can burn up to 3 times more calories than any traditional exercise bike.

#### Usability

Your client can easily choose specific programs or set the parameters for each session based on his needs.

#### Tracking


The client can easily track his performance in real time.

#### Effect

BodyTone Lab Horizontal takes fitness to a new level! Not only maximizes the efficiency of the workouts, but also rejuvenates the body with a long-lasting efficacy. Features Outperforming expectations at first sight

A stylistic modern design of a dynamic capsule bike that combines innovative technical solutions:

- Touch screen display
- Luxury leather
- Under-pressure vacuum massage
- Anti-gravity exercises
- Infrared heating
- Collagen lamp
- Color and aroma therapies
- A multivision virtual with simulator and mp3 music



The patient can enjoy a futuristic exercise experience on a luxury bike chair in a comfortable lying down position



### Product Details and Specifications

Infrared	<ul> <li>Temperature adjustment: 25° C - 50° C</li> <li>Two infrared radiators: 350W each</li> <li>Protection against scalding</li> <li>The range of infrared radiation: 1.4 μm - 14μm (close to the solar radiation)</li> </ul>
Vacuum	5 predefined programs: 0 - 30mBar A function to adjust the vacuum manually A vacuum pump - 1500W
Chromotherapy	Based on LED technology Wireless remote control allowing to choose the colour Predefined programs for automatic colour change
Aromatherapy	An automatic dosing station with essential oils for the THERMAL operation mode An option to turn on/off the aromatherapy during the training
Air ionization	Negative ion emission: 9,000,000 ions/cm3
Ozone	Ozone emission: 0.05ppm
Collagen	A single LED lamp - 30W Wavelength - 620nm An option to turn on/off during the training



## BODYTONE



LASERVILLE Medical Equipment Trading LLC

+971 55 386 5007

info@laserville.net

 $\square$ 

www.laserville.net